

Made in France



HORIZON[®]
P5

Portable Oxygen Concentrator

Horizon[®] proposes premium and highly reliable range of oxygen concentrators designed for oxygen therapy patients who want to maintain their lifestyle.

Bluetooth connected

Light and intuitive

High O₂ purity

SCALEO
MEDICAL

BENEFITS

LIGHTWEIGHT
& COMPACT



20 SETTINGS FOR
PRECISE DOSAGE

from 0,25 to 5

EASY
MAINTENANCE



HORIZON[®]
P5



3 YEARS DEVICE'S
WARRANTY**

BLUETOOTH
CONNECTED*



LONG LIFE
BATTERY

* IF ACTIVATED ** BATTERY, CHARGER AND SIEVE BEDS WARRANTED FOR 1 YEAR

HORIZON[®]
P5

pack includes:

- Horizon[®] P5 portable concentrator
- Single battery or double battery
- Carry bag
- AC/DC power supply 100-240 V, 50-60 Hz
- User Manual *not shown in picture*



Intensive use



Designed & manufactured in France

ACCESSIBILITY



Simple & intuitive



High oxygen purity



Easy to access filters



Charging and maintenance access

TECHNICAL FEATURES

Dimensions w x d x h	216 x 88 x 220 mm*
Weight	2,54 kg*
Sound level	44 dB at setting 2
Oxygen purity	> 87 %
Pulse flow in mL/min	Up to 1050 mL/min
Oxygen flow rate setting	20 flow settings (from 0,25 to 5)
Battery life	Up to 5,5 hours on setting 1**
Recharge time	Up to 2,5* hours and 4,5** hours

Power consumption	85 W
Power compatibility	100 - 240 VAC / 50 HZ 60 HZ
Environmental tolerance	Temperature +5 to +40°C Relative humidity : 95 %
Storage environmental tolerance	Temperature: -25 to +70°C Relative humidity: 0 % - 95 %

* with simple battery
** with double battery

SCALEO
M E D I C A L

Designed and manufactured by
SCALEO Medical - 107 rue Joe Dassin, Parc 2000 - 34 080 Montpellier - France
Tel. +33 467 72 94 86 - Fax +33 499 74 01 79
www.scaleomedical.com - info@scaleomedical.com

CE
0459